

## ***In-situ* observation in aberration corrected TEM**

T. Oikawa<sup>1)</sup> and I. Ohnishi<sup>2)</sup>



1) JEOL Asia Pte Ltd., #01-12 Corporation Place, 2 Corporation Road 618494, Singapore

2) JEOL Ltd., 1-2 Musashino 3-Chome, Akishima, Tokyo 196-8558, Japan

Corr. Author (e-mail): oikawa@jeolasia.com.sg

*Keywords:* *In-situ*, aberration correction, cold FEG, catalyst, electro plating

### **Abstract**

*A various experimental results in In-Situ observations which were combined with aberration corrected transmission electron microscopy (Cs-corrected TEM) are introduced. In this lecture, the features and the performances of newly developed Cs-corrected TEM and a various kind of In-Situ specimen holders are introduced in detail.*

### **Introduction**

Owing to the aberration corrected TEM/STEM, atomic resolution observation and analysis are commonly used, today [1]; however, the TEM specimen is kept under high-vacuum condition. It is very important to investigate the structures analysis of materials under practical operating environment; therefore the *In-Situ* observation technology is becoming very important. Recently, a various kinds of *In-Situ* specimen holders have been developed; the *In-Situ* observation can be introduced into the market. And also, JEOL developed Cs-corrected TEM with wide gap objective polepiece for accepting the *In-Situ* holders and offering higher efficiency X-ray signal collection.

### **Results and discussions**

Figure 1 shows an appearance of newly developed Cs-corrected TEM (*JEM- ARM300FII*). The cold field emission gun (cold-FEG) is employed and owing to a small energy spread and a high-coherency of electrons of cold-FEG, a TEM lattice resolution of 60 pm (Cs-corrected with WGP at 300 kV) is performed. And dual type of largest size of EDS detector (158 mm<sup>2</sup>, each) can be installed into the wide-gap objective polepiece (WGP); thus, the X-ray collection solid angle becomes very large as 2.21 sr. for analytical capability.

Figure 2 shows an example of *In-Situ* observation of *Ir/CeO<sub>2</sub>* catalyst redox reaction obtained at specimen temperature of 600 deg. C. with gas environment holder (*Protochips Atmosphere*) and *ARM300FII*, and gases were transferred as  $H_2 - O_2 - H_2$ . The *Ir/CeO<sub>2</sub>* catalyst redox reaction was observed at atomic resolution.

Figure 3 shows another example of *In-Situ* observation of electrochemical reaction (copper electro plating) obtained with liquid environment holder (*Protochips Poseidon*) and *ARM300FII*. The copper plating layer growing in a cupric sulfate solution is observed dynamically. Owing to *In-Situ* specimen holders and aberration corrected TEM combination, the *In-Situ* observation is realized.

### **References**

[1] C. Ricolleau et al.; *Microscopy* 62 (2013) 283 – 293.




Fig. 1. Appearance of newly developed Cs-corrected TEM (*JEM-ARM300FII*).

Fig. 2. *Ir/CeO<sub>2</sub>* catalyst redox reaction @600 deg C. The gases were transferred as *H<sub>2</sub> – O<sub>2</sub> – H<sub>2</sub>* in gas environment holder.




Fig. 3. Electrochemical reaction (copper electro plating) The copper plating layer growing in a cupric sulfate solution is observed dynamically with liquid environment holder (*Protochips Poseidon*) and *JEM-ARM300FII*. Detector: *Gatan OneView Insitu*.