# INDIAN INSTITUTE OF TECHNOLOGY DEPARTMENT OF CIVIL ENGINEERING

### **Enquiry letter for purchase of SMPS**

Sub: Quotation for supply, installation, commissioning and training of Scanning Mobility Particle Sizer (SMPS)

Reference: IITK/CE/2015/1029 Dated: October 1, 2015

Sir / Madam.

With reference to the subject mentioned above, you are invited to submit the quotation in a sealed cover. Configuration/Specification are given below:

# **Specifications of Scanning Mobility Particle Sizer:**

**Principle of operation:** The Scanning Mobility Particle Sizer is based on the principle of the mobility of a charged particle in an electric field. Particles entering the system are neutralized (through a bipolar diffusion charger) such that they have a Fuchs equilibrium charge distribution. They then enter a Differential Mobility Analyser (DMA) where the aerosol is classified according to electrical mobility, with only particles of a narrow range of mobility exiting through the output slit. The particle flow is injected at the outside edge of the DMA, particles with a positive charge move across the sheath flow towards the central rod, at a rate determined by their electrical mobility. This monodisperse distribution then goes to a Condensation Particle Counter which determines the particle concentration at that size.

SMPS Settings and Requirements High resolution data  $64 \text{ channels/decade, total } \geq 165$ Broad size range from 2.5 nm to  $\geq 1,000 \text{ nm}$ Wide concentration range minimum up to  $10^7 \text{ particles/cm}^3$ 

Data averaging (Scans per Sample) 1 to 999, user-selectable

Aerosol Flow rate  $0.2 \text{ to } \ge 3 \text{ l/min, user-adjustable}$  Sheath Flow rate 2 to 30 l/min, user-adjustable

Scan time < 10 sec

Working Fluid alcohol

Operating temperature 10 to 35°c

Storage temperature 0 to 40°c

Aerosol-inlet temperature 10 to 35°c

Humidity 0 to 90%, noncondensing

Pressure 75 to 105 kPa

Operation system should be capable of using

with computer and without computer

Aerosol neutralizer Bi-polar diffusion charger/

Radioactive/X-Ray

# **Data Logging**

Via attached Pc running Microsoft® Windows® (Windows® 7 32/64 bit compatible) as well as inbuilt data storage should be possible

## **Display**

LCD display on DMA controller and CPC

#### **Communications**

RS-232, USB, and Ethernet for full control

#### **Inlet accessory**

Single-stage, inertial impactors (choice of three impactors, each with a different cut size)

## **Power Requirements**

Instrument should be capable of working at Indian voltage standards.

# Support/Service capability

Supplier should have supplied min. 5 systems in India and provide details of their Engineers capable of after sales Service support.

## The quotation should have the following details:

- 1. Indicate item-wise pricing on FOB/CIF basis
- 2. Technical specifications in detail
- 3. Technical bid and Price bid to be sealed separately (Two-bid system)
- 4. Warranty period
- 5. User list of similar system supplied in India.
- 6. Maximum educational discount considering end use for research and teaching
- 7. Payment terms
- 8. Proprietary Certificate, if applicable
- 9. Support/Service capability in INDIA
- 10. Comprehensive AMC prices should be quoted separately
- 11. In case of CIP/CIF price, insurance from your w/h to our w/h in IIT Kanpur must be arranged by the bidder/supplier
- 12. Any other relevant details

#### Terms and condition:-

- 1. Sealed Quotation must reach the undersigned on or **before October 10, 2015** (Extended up to October 30, 2015).
- 2. Prices should be in USD and CIF Delhi.
- 3. Our Institute is partially exempted from custom duty.
- 4. The final selection will be made based on weights given to technical merit and pricing as 70% and 30% each, respectively.

Dr. Sachchida Nand Tripathi Professor Department of Civil Engineering Indian Institute of Technology-Kanpur Kanpur – 208016

Tel. 0512 – 2597845 Fax No. 0512-2597395 E-mail: snt@iitk.ac.in